**EN**   READ CAREFULLY BEFORE USE

### Light-cured Dental Restorative Material

BEAUTIFIL II Gingiva is a light-cured radiopaque restorative material, designed to restore cervical defects as well as to cover receded gum.

**SHADES**

G-DP (G-Dark Pink), G-LP (G-Light Pink), G-Br (G-Brown), G-Or (G-Orange), G-V (G-Violet)

**COMPOSITION**

Bis-GMA, TEGDMA, S-PRG filler (fluoroboroaluminosilicate glass), Polymerization initiator, Pigments and others

**INDICATIONS**

➢ Class V restorations (Cervical caries, root caries, wedge-shaped defects)
➢ Covering of exposed or discolored necks of teeth
➢ Restoration of esthetic appearance of gingival areas

**PRECAUTIONS**

- Do not use this product on patients with known allergies to this material and/or methacrylate monomer.
- Operators with known allergies to this material and/or methacrylate monomer must not use this product.
- In case of inflammation or other allergic reactions in either patient or operator while using this product, immediately discontinue use and seek medical advice.
- Use medical (dental) gloves and protective glasses to avoid sensitization to this product. In case of accidental contact with oral soft tissue or skin, immediately blot with alcohol moistened cotton ball, and rinse with plenty of water. Medical (dental) gloves provide some protection against direct contact with this material; however some monomers can penetrate through the gloves very quickly to reach the skin. Therefore, should this material contact even the gloves, immediately remove them and wash hands with plenty of water. In case of contact with eyes, flush immediately with plenty of water and seek medical advice.
- This product is intended for use by dental professionals only.

**DIRECTIONS FOR USE**

- Cleaning of tooth surface
  - Thoroughly clean the tooth surface to remove plaque/surface stains.
- Shade selection
  - Select appropriate shade using BEAUTIFIL II Gingiva shade guide.
- Cavity preparation
  - Prepare the cavity following the conventional method. In case of non-carious defects, roughen the tooth surface to enhance mechanical and chemical bond. Use a gingival retraction cord where necessary.
- Isolation
  - Use a rubber dam etc to isolate the restoration site.
- Pulp protection
  - In the case of pulp exposure or deep cavities with close proximity to pulp, cap with calcium hydroxide etc.
- Adhesive system
  - Apply dental adhesive following the manufacturer’s instructions for use, such as SHOFU’s BeautiBond, FL-BOND II.
- Restoration
  - Dispense adequate amount of BEAUTIFIL II Gingiva from the syringe onto a paper pad. Place the dispensed material into the cavity or tooth surface and contour to desired shape with a suitable instrument. For deeper application, place material in 1.5 mm or less increments.
    - When storing these materials in a refrigerator, allow them to reach room temperature before use.
    - Hold the syringe securely to remove the cap.
    - The paste dispensed onto the paper pad should be shielded from ambient light with light-shielding cover to avoid premature polymerization of the material.
    - Tightly cap the syringe after each use.
  - Light-curing
    - Light-cure using a standard dental light-curing unit.
    - In deeper application, light-cure each layer.
    - In cases where the restored area is large, light-cure in segments.
    - Specifications of light-curing units recommended for this material and their standard curing time

Light-curing unit	Halogen	LED
Light-curing time	40 sec.	20 sec.
Wavelength	400-500 nm	440-490 nm
Light intensity	≥ 500 mW/cm²	≥ 1000 mW/cm²

➢ Use a dental light-curing unit with specifications equivalent to those described above. When using the light-curing units with performance scale other than these recommendations, follow the instructions of each unit to light-cure BEAUTIFIL II Gingiva properly.
➢ Operators must use protective glasses or light shielding plate to avoid direct exposure to the curing light. Operators must also protect patients' eyes from the curing light.
➢ Disinfect the light-curing unit with alcohol after each use to avoid cross contamination.
9. Finishing and polishing
After light-curing, contour, finish and polish the restoration with OneGloss or Super-Snap finishing and polishing system.

**NOTE ON USE**

- Do not use eugenol containing materials since they might inhibit polymerization or adhesion of the material.
- Follow the instructions for use of each dental material, instrument or equipment to be used in conjunction with this product.
- Do not use this product for any purposes other than specifically outlined in these instructions for use.
- Use this product within the expiration date indicated on the package and container.
(Example ☞ YYYY-MM-DD → Year-Month-Date of the expiration date)

**STORAGE**

Store at room temperature (1-30 °C/34-86 °F). Avoid high temperatures and high humidity. Keep away from direct sunlight and any source of ignition.

***CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.***

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